

## ChIP-seq with *N. furzeri* embryos

### Materials

Proteinase K solution (Invitrogen 25530-049)

RNase A (Sigma R6513, dissolved in 10mM Tris-HCL, pH7.5, 15mM NaCL at concentration of 10 mg/mL)

Complete™, Mini, EDTA-free Protease Inhibitor Cocktail (Sigma 11836170001)

Dynabeads™ M-280 Sheep Anti-Rabbit IgG (Thermo Fisher 11203D)

Paraformaldehyde, 16% w/v aq. soln., methanol free, 10 x 10mL (Alfa Aesar 43368-9M)

Rabbit anti-H3K27me3 antibody (Millipore, 07-449 or Active motif 39155)

Rabbit anti-H3K4me3 antibody (Active Motif 39159)

DTT 1M (Thermo Fisher Scientific)

GlycoBlue 15 mg/mL (Ambion, AM9515)

16% Formaldehyde solution (w/v), Methanol-free (Thermo Scientific, 28906)

### Recipes

*All buffer should be filter-sterilized and store at 4°C.*

| Buffer A2 (250ml, w/o PI) | Reagent   | Volume  | Final  |
|---------------------------|---|---------|--------|
|                           | 1M HEPES pH 7.5   | 3.75 ml | 15 mM  |
|                           | 5 M NaCl  | 7 ml    | 140 mM |
|                           | 0.5M EDTA   | 0.5 ml  | 1 mM   |
|                           | 0.5M EGTA   | 250 ul  | 0.5 mM |
|                           | Triton X100   | 2.5 ml  | 1%     |
|                           | 10% sodium deoxycholate                                     | 2.5 ml  | 0.1%   |
|                           | 10% SDS   | 25 ml   | 1%     |
|                           | 10% N-lauroylsarcosine                                      | 12.5 ml | 0.5%   |
|                           | (Add DTT (final 0.5mM) and proteinase inhibitor before use) |         |        |

| Low-salt buffer (500 ml) | Reagent            | Volume | Final  |
|--------------------------|--------------------|--------|--------|
|                          | 1M Tris-HCl pH 8.0 | 10 ml  | 20 mM  |
|                          | 5M NaCl            | 15 ml  | 150 mM |
|                          | 0.5M EDTA          | 2 ml   | 2 mM   |
|                          | 10% Triton X-100   | 50 ml  | 1%     |

| High-salt buffer (500 ml) | Reagent            | Volume | Final  |
|---------------------------|--------------------|--------|--------|
|                           | 1M Tris-HCl pH 8.0 | 10 ml  | 20 mM  |
|                           | 5M NaCl            | 50 ml  | 500 mM |
|                           | 0.5M EDTA          | 2 ml   | 2 mM   |
|                           | 10% Triton X-100   | 50 ml  | 1%     |

| LiCl buffer (500 ml) | Reagent                 | Volume | Final  |
|----------------------|-------------------------|--------|--------|
|                      | 1M Tris-HCl pH 8.0      | 10 ml  | 20 mM  |
|                      | 5M LiCl                 | 25 ml  | 250 mM |
|                      | 0.5M EDTA               | 1 ml   | 1 mM   |
|                      | 10% Triton X-100        | 50 ml  | 1%     |
|                      | 10% NP-40               | 5 ml   | 0.1%   |
|                      | 10% sodium deoxycholate | 25 ml  | 0.5%   |

### Harvest cells [Day 1]

Wash 20-40 fin blastema with PBS twice @RT. Cut samples into small pieces on ice for 5 min in a well of 24 well plate.

Add 1 ml 1% FA in PBS and crosslink for 12 min. Quickly transfer crosslinked samples into 1.5 ml tubes and centrifuge for 1 min at 10000 rpm to quickly collect cells.

Remove fixation buffer and add 2.5M glycine to 0.125M final to quench crosslink.

Wash with ice-cold PBS/0.125M glycine twice and once with Buffer A2.

*You can keep the pellet at -80°C after a snap freeze.*

### Prepare beads [Day 2, just before shear chromatin]

(Vortex the bottle of beads well before use, and below using ab4729 as an example)

Add 130µl Rabbit IgG Dynal beads to a 1.5ml microfuge safe-lock tube for each ChIP.

Wash beads with PBS twice.

Add 700µl PBS-0.5% BSA with 10µg ab4729 antibody, and incubate 2-6h @ 4°C.

**(Please check the affinity of beads with the antibodies you use. For Protein A/G, 50µl beads for 10µg antibody)**

### Shear chromatin [Day 1-2]

Resuspend cells ( $\approx 0.5-1 \times 10^7$ ) with 700-800µl Buffer A2 per tube.

Sonicate 8-10  $\times$  10s @ power 3 with Microson ultrasonic cell disrupter till the size of DNA from 150 to 1000 bp. (at power setting #5, 30 seconds sonication followed by 30 seconds interval, 15 rounds total.

Centrifuge 10 min @ max speed @ 4°C and, keep supernatant for ChIP.

### ChIP [Day 2]

Wash beads with PBS-0.5% BSA twice.

Add 700µl sheared chromatin, and incubate 6h or O/N @ 4°C. **(Leave 30ul for input)**

Wash beads-chromatin twice with low salt buffer, twice with the high-salt buffer, and once with LiCl buffer @ cold room. (Keep the supernatant to check sonication efficiency)

### **Purify DNA [Day 3]**

Add 150µl elution buffer, 150µl TE, and 4µl RNase, incubate at 37°C, with manual shaking every 10 minutes for 30 min total (For input, add 120µl elution buffer, 150µl TE, and 4µl RNase)

Add 2µl Proteinase K and incubate at 65°C with manual shaking every 20 minutes for 6h total.

Extract with 300µl phenol-ch-iso for (25:24:1), vortex 1 min then centrifuge 5 min at 20000g.

Take out 280-300µl upper phase and transfer to a new microfuge tube containing 12µl of 5M NaCl (final = 200mM) and 6µl glycogen (final = 30µg)

Mix well and add 750µl cold (4°C) 100% EtOH. Keep @ -80°C for 1h or @-20°C O/N then centrifuge at max speed for 30min at 4°C and wash once with 1ml 75% EtOH.

Air-dry 30min at RT, and then resuspend pellets in 40µl ddH<sub>2</sub>O for 1h @ 4°C

Check sonication condition by running a 1.7% DNA gel.

### **ChIP elution buffer (50ml)**

10 mM Tris, pH8.0

1 mM EDTA

1% SDS

250 mM NaCl

v1.0: Chi-Kuo Hu, 2021/09/01

Adapted from the working protocol of

Wei Wang @ Alejandro Sánchez Alvarado's Lab, Stowers, 2019/02/21