

## *in situ* HCR v3.0 with *N. furzeri* embryos

### Materials

- RNase free 8% Paraformaldehyde (PFA) (Thermo Fisher Scientific, 47347)
- 10× Ca<sup>2+</sup> and Mg<sup>2+</sup> free phosphate-buffered saline (PBS) (Research Products International, P32200)
- Tween 20 (Fisher Chemical, BP337)
- Proteinase K, RNA grade, 20 mg/mL (Invitrogen, 2553004)
- 20× Sodium chloride sodium citrate (SSC) (Sigma-Aldrich, S6639)
- Methanol (Fisher Chemical, A411)
- Glycerol (Fisher Scientific, BP229)
- Hoechst (Thermo Fisher Scientific, H3570)

### Recipes

<b>PBS (50 mL)</b>	<b>Reagent</b>	<b>Amount</b>	<b>Final</b>
	10× PBS	5 mL	1× PBS
	Fill up to 50 mL with ultrapure H <sub>2</sub> O		
<b>PBST (50 mL)</b>	<b>Reagent</b>	<b>Amount</b>	<b>Final</b>
	10× PBS	5 mL	1× PBS
	10% Tween 20	500 µL	0.1% Tween 20
	Fill up to 50 mL with ultrapure H <sub>2</sub> O		
<b>4% PFA (20 mL)</b>	<b>Reagent</b>	<b>Amount</b>	<b>Final</b>
	8% PFA	10 mL	4% PFA
	10× PBS	2 mL	1× PBS
	ddH <sub>2</sub> O	8 mL	
<b>Proteinase K solution (1 mL)</b>	<b>Reagent</b>	<b>Amount</b>	<b>Final</b>
	20 mg/mL proteinase K	1.5 µL	30 µg/mL
	Fill up to 1 mL with PBST		
<b>5× SSCT (40 mL)</b>	<b>Reagent</b>	<b>Amount</b>	<b>Final</b>
	20× SSC	10 mL	5× SSCT
	10% Tween 20	400 µL	0.1% Tween 20
	Fill up to 40 mL with ultrapure H <sub>2</sub> O		

### The storage conditions for HCR kits

- Store HCR probe sets, HCR amplifiers, HCR probe hybridization buffer, and HCR probe wash buffer at -20 °C.
- Store HCR amplification buffer at 4 °C. Please mix well and make sure there are no precipitates in the buffer when you bring it to room temperature before use.
- On the bench top, keep stock solutions on ice.
- Make sure all solutions are well mixed before use.

### Methods

**This protocol is carried out at room temperature (RT) unless otherwise indicated.**

#### Fixation [Day 1]

1. Collect embryos and incubate at 27 °C in a petri dish with Ringer's solution.
2. Transfer the embryos into a clean glass dish with PBS. Use clean fine tip forceps to dechorionize the embryos and remove the enveloping layer membrane.
3. Transfer 5-10 embryos to a 1.5 mL Eppendorf tube and remove excess PBS using the 200 µL pipette.
  - ***Note:** Be aware that embryos might stick and stay in the pipette during transfer. Be careful not to accidentally discard them with the pipette tip. Early-stage embryos in the tube can be briefly spun down in a microcentrifuge for 3 s before any solution changes.*
4. Fix embryos in 250 µL 4% PFA made in PBS. Incubate for 2-3 h at RT or overnight at 4 °C.
  - ***Caution:** Use PFA with extreme care, as it is a hazardous material.*
  - ***Note:** Use fresh PFA and cool to 4 °C before use to avoid increased autofluorescence.*
5. Wash embryos in 250 µL PBS 3 times for 5 min each to stop the fixation.
  - ***Note:** Avoid using calcium chloride and magnesium chloride in PBS, as this leads to increased autofluorescence in the embryos.*

6. Dehydrate and permeabilize embryos in a series of graded 250  $\mu$ L MeOH/PBS washes sequentially for 5 min each:
  - (a) 25% MeOH / 75% PBS (e.g. 10 mL MeOH and 30 mL PBS for 40 mL solution)
  - (b) 50% MeOH / 50% PBS (e.g. 20 mL MeOH and 20 mL PBS for 40 mL solution)
  - (c) 75% MeOH / 25% PBS (e.g. 30 mL MeOH and 10 mL PBS for 40 mL solution)
  - (d) 100% MeOH.
7. Store at -20°C overnight before use.

**Pause Point:** *Samples in 100% MeOH can be stored for six months at -20 °C.*

### **Prepare to stain [Day 2]**

8. Rehydrate embryos in a series of graded 250  $\mu$ L MeOH/PBST washes sequentially for 5 min each:
  - (a) 75% MeOH / 25% PBST
  - (b) 50% MeOH / 50% PBST
  - (c) 25% MeOH / 75% PBST
  - (d) 100% PBST.
9. Wash embryos in 250  $\mu$ L PBST 2 times for 5 min each.
10. Treat embryos in 250  $\mu$ L proteinase K. Wash embryos in 250  $\mu$ L PBST 2 times without incubation. Postfix embryos with 250  $\mu$ L 4% PFA for 20 min. Wash embryos in 250  $\mu$ L PBST 5 times for 5 min each.
  - **Note:** *Skip the proteinase K treatment and the following postfixation step for embryos before pharyngula stage, including embryos in diapause.*
  - **Note:** *Proteinase K concentration and treatment time should be reoptimized for each batch of proteinase K or for samples at a different developmental stage.*

**Pause Point:** *It is reported that samples in PBST can be stored for up to 4 days at 4 °C. However, the hybridization step is recommended to be started immediately after staining preparation.*

### **Detection Stage [Day 2]**

11. For each sample, transfer 10 embryos to a 1.5 mL tube.
12. Pre-hybridize embryos in 250  $\mu$ L probe hybridization buffer at 37 °C metal bath for 30 min.

- **Caution:** Probe hybridization buffer contains formamide, a hazardous material.
  - **Note:** Pre-heat enough probe hybridization buffer (250  $\mu$ L for pre-hybridization + 250  $\mu$ L for hybridization per tube) to 37°C before use.
  - **Note:** Early-stage embryos tend to be suspended in the probe hybridization buffer. Briefly spinning embryos down will help the solution change.
13. Prepare probe solution by adding 1  $\mu$ L of 1  $\mu$ M probe mixture to 250  $\mu$ L of probe hybridization buffer at 37 °C. For multiplexing, add 1  $\mu$ L of each 1  $\mu$ M probe mixture to the same 250  $\mu$ L of probe hybridization buffer at 37 °C.
- **Note:** For a 1  $\mu$ M concentration stock, DNA oPools (50 pmol/oligo) are dissolved in 50  $\mu$ L TE buffer and stored at -20 °C.
  - **Note:** The standard reaction concentration of each HCR probe is 4 nM. The HCR probe set scale depends on the reaction volume and the number of experiments. For example, for a 250  $\mu$ L reaction volume, each experiment will require 1 pmol of each probe, so a 50 pmol probe set scale corresponds to 50 experiments.
14. Remove the pre-hybridization solution and add 250  $\mu$ L probe solution to each tube.
15. Incubate embryos overnight (12–16 h) at 37 °C.

### Detection Stage [Day 3]

16. Remove excess probes by washing embryos in 250  $\mu$ L probe wash buffer at 37 °C 4 times for 15 min each.
- **Caution:** Probe wash buffer contains formamide, a hazardous material.
  - **Note:** Pre-heat enough probe wash buffer (1000  $\mu$ L per tube) to 37°C before use.
  - **Note:** Early-stage embryos tend to suspend in the probe wash buffer. Brief centrifugation can make solution change easier.
17. Wash embryos in 250  $\mu$ L 5 $\times$  SSCT 3 times for 5 min each.

**Pause Point:** Samples can be stored in SSCT for hours. However, the amplification step is recommended to be started immediately after the hybridization step.

### Amplification Stage [Day 3]

18. Pre-amplify embryos in 250  $\mu$ L amplification buffer for 30 min.
- **Note:** Equilibrate enough amplification buffer (250  $\mu$ L for pre-amplification + 250  $\mu$ L for hairpin solution per tube) to RT before use.

19. Add 5  $\mu$ L hairpin h1 and 5  $\mu$ L hairpin h2 from the 3  $\mu$ M stock in separate 0.2 mL PCR tubes. For multiplexing, add 5  $\mu$ L individual hairpin h1 and 5  $\mu$ L individual hairpin h2 from the 3  $\mu$ M stock in separate 0.2 mL PCR tubes. Snap cooling the hairpin tubes in a PCR machine as follows:
  - $\left\{ \begin{array}{ll} \text{(a) } 95^{\circ}\text{C} & 1\text{min } 30\text{s} \\ \text{(b) } 0^{\circ}\text{C} & 5\text{min} \\ \text{(c) } 26^{\circ}\text{C} & 30\text{min} \end{array} \right.$
- **Note:** Hairpins should be protected from light as much as possible.
- **Note:** HCR hairpins h1 and h2 are provided in hairpin storage buffer ready for snap cooling. All hairpins should be snap cooled in separate tubes.
- **Note:** The standard reaction concentration of each HCR amplifier hairpin is 60 nM. The HCR amplifier scale depends on the reaction volume and the number of experiments. For example, for a 250  $\mu$ L reaction volume, each experiment will require 15 pmol of each hairpin, so 1 nmol amplifier scale corresponds to 66 experiments.
20. Prepare hairpin solution by adding snap-cooled h1 hairpins and snap-cooled h2 hairpins to 250  $\mu$ L amplification buffer.
21. Remove the pre-amplification solution and add 250  $\mu$ L the hairpin solution to each sample.
  - **Note:** Keep in dark for the rest of the protocol.
22. Incubate the embryos overnight (12–16 h) in a box in a dark drawer.

#### **Amplification Stage [Day 4]**

23. Remove excess hairpins by washing embryos in 250  $\mu$ L 5 $\times$  SSCT 3 times for 5 min each.

**Pause Point:** Samples in SSCT can be stored at 4  $^{\circ}$ C for a week or longer protected from light before microscopy.

#### **Microscopy [Day 4]**

24. Remove 5 $\times$  SSCT and add 250  $\mu$ L of 5  $\mu$ g/mL Hoechst staining solution for around 3 h.
25. Clear embryos by submerging them in a series of 500  $\mu$ L glycerol/SSCT solutions sequentially for 5 min each:
  - (a) 30% glycerol / 70% SSCT

- (b) 50% glycerol / 50% SSCT
- (c) 70% glycerol / 30% SSCT

**Pause Point:** *Samples in 70% glycerol / 30% SSCTPBST can be stored at 4 °C for 3 months or longer protected from light before microscopy.*

26. Place embryos on a clean glass slide without excess glycerol solution. Gently stage the embryos for imaging. Apply a small amount of vacuum grease to each corner of a clean coverslip (not over the edges) and place the coverslip on top of the embryos. Gently tap the corners to seal.
  - **Note:** *Embryos might deform when over-pressing the coverslip. Make sure the coverslip is horizontal or else the embryo will drift after mounting.*
27. Acquire images using a confocal microscope.

Adapted from the working protocol of

Benjamin Martin's Lab, Stony Brook University, 2021/03/02

v1.0: Lu Jia, 2022/08/05